







## KB MEDICAL(GROUP), INC.

Our Headquarter

Add.: 9650 Research Drive, Suite A, Irvine, California 92618, USA

Tel: 1-949-333-3785

Email: info@kbmed.com

Website: www.kbmed.com

Our production unit

Jiangsu Kangbao Medical Equipment Co. Ltd.

Add.: No. 78, North Suzhong Road, Baoying, Yangzhou city, Jiangsu, China 225800

Our commercial office

Hangzhou Best Valley Medical Technology Co. Ltd.

Add.: Room 1105, Tower A, Wufutianxinglong Mansion,

No. 358 East Fengqi Road, Jianggan district, Hangzhou, China 310020

Tel: 0086-571-87987813, 87975379, 87983593

Email: susan@bestvalleymedical.com

Website: www.kbmed.com

















## **KB Medical Group**

KB Medical Group was founded in Irvine, California.

A professional company specialized in producing innovative Safety Medical Products to protect health care workers against sharps injury. KB Medical Group is committed to address these issues by providing the most effective safety-engineered therapy devices for optimum safety in the hospitals.

In the past several years, KB Medical Group focus on working on a wide variety of innovative, high-quality enteral feeding products and accessories uniquely designed for delivering nutrition and medication to pediatric and adult patients. This new standard ENFit® connection has been introduced in-line with ISO80369-3 has a unique enteral-specific design that help reduce the risk of enteral tube feeding misconnections and improve the patient safety.

#### **KB Medical Group in the world**


- Our Headquarter: Irvine, California, USA.
- Our production unit: Jiangsu Kangbao Medical Equipment Co. Ltd. is our professional manufacturer in Jiangsu province, China (1 hour away Shanghai) for more than 30 years.
- Our commercial office: Hangzhou Best Valley Medical Technology Co. Ltd. is our Commercial Center in China which is located at Hangzhou, Zhejiang province (40minutes away Shanghai), the center in charge of the business communication with customers, arrange the customers for visiting plan, QA audit, also provide service for the registration documents and documentary works.



#### **Key facts of our KB Medical Group**



#### **Quality, Certificate, Innovation**



### Key facts of our production unit



40 sets Injection molding machines

16 sets Auto printing lines for syringes

8 sets Auto-assembling machines for syringes

8 sets Blister packing machines

8 sets
Auto-loading packing machines for syringe (1/3/5/10/20/60/100ml)

1 set Auto-assembling machine for infusion set

1 set Auto-assembling machine for scalp vein set

3 sets

Auto-assembling machines for hypodermic needle

3 sets

Modernization testing lab for Physical, Chemical and Biological

3 sets EO Sterilization chamber





► Hypodermic needle workshop



#### **CATALOGUE CATEGORIES**



#### Infusion therapy



Wide range of a safety way for injection of vein for medication infusion. IV cannula

- Safety scalp vein set
- Scalp vein set Infusion set
- IV burette
- Safety IV cannula IV cannula dressing
- Sterile IV start kit
- Extension tube
  - Three way stopcock

#### Page 1-9

#### Syringe and needle

A complete range of syringes for injection or harvesting of biological material with safety protection.

- Disposable syringe with needle
- Disposable syringe with safety needle
- Safety syringe
- Safety needle
- Disposable hypodermic needle
- Auto-Disable syringe for fixed-dose immunization
- Pre-filled flush syringe
- Dental needle
- Irrigation needle
- Suture removal Set
- Veterinary needle
- Hard pack syringe



Page 10-21

#### Hazardous drug safety



Various drug preparation needles for safety drug admixture and disposal of used sharps.

- Aspiration needle with safety shield Blunt needle with/without filter
- Curved needle
- Side hole needle
- Safety box for the disposal of used sharps
- Needle free system
- Sharp container
- Endoscope transport pad
- Ultimate sharps safety stations (needle counter)
- Safety medical products for NEEDLE STICK PREVENTION

Page 22-28

#### Safety medical products

Wide range of safety medical products to protect health care workers against sharps injuries.

- Safety scalp vein set
- Safety IV cannula
- Safety needle Safety syringe
- Auto-Disable syringe
- Safety blood lancet
- Safety insulin pen needle
- Safety insulin syringe
- Safety fistula needle
- Safety blood collection needle
   Safety drug preparation needle
  - Ultimate sharps safety stations
  - Needle free system
  - Safety box for the disposal of used sharps
  - Sharp container



Page 28

#### Nutrition therapy



Various components for nutrition therapy which is compliant with the ISO80369-3 to help reduce the risk of misconnections.

- ENFit®syringe
- The ENFit® connector
- ENFit® Feeding set
- Oral syringe
- Irrigation syringe

#### Specimen collection



Page 33-36

Blood collection needle Safety blood collection needle

- Our safety blood collection system help eliminate accidental needle stick injuries when collection blood sample.
  - Safety lancet
  - Blood collection tube

#### Anesthesia & pain management

Wide range of products for subdural, lower thorax and lumbar spinal puncture.

- Spinal needle
- High pressure manifold
- Anesthesia kit
- NRFit™ syringe



Page 37-39

#### Transfusion therapy



Page 40-42

Complete range of safety way to collect, preserve, filtration and administration for human blood and blood components.

- Gravity blood administration sets
- Blood bag
- Cord blood bag

#### Urinary care

We provide solutions to help improve patient safety and infection prevention in urinary care.

- Urine drainage bag
- Leg bag
- Urine meter
- Irrigation tray Irrigation syringe



#### Renal therapy



Innovative, high performance products and services for use in dialysis, provide a basis which targets improved clinical outcomes for the patient, while minimizing costs and risks for the healthcare providers.

- Fistula needle
- Blood line
- Safety fistula needle
   Hemodialysis care kit







Safety insulin syringe

Page 47-48

- Insulin syringe
- Insulin pen needle
- Safety insulin pen needle • Sharp disposal container for insulin pen needle
- Alcohol swab



Page 49-52

### Infusion therapy



### Safety scalp vein set











#### Features

- ► The ergonomic design is ideal for clinician handling and patient comfort.
- Offering an easy grip.
   The design to minimize
- The design to minimiz slips and touch contamination.









Co	de	Gauge	Size and length		Color
U	V			V V	U
SSVS27/	5VS27	27G×5/8"-3/4"	0.4×16-19mm		Grey
SSVS26/	5VS26	26G×5/8"-3/4"	0.45×16-19mm		Brown
SSVS25/	5VS25	25G×5/8"-4/5"	0.5×16-20mm		Orange
SSVS24/	5VS24	24G×5/8"-4/5"	0.55×16-20mm		Purple
SSVS23/	5VS23	23G×3/4"-1"	0.6×19-25mm		Blue
SSVS22/	5VS22	22G×1"	0.7×25mm		Black
SSVS21/	5VS21	21G×1-1/8"	0.8×28mm		Green
SSVS20/	5VS20	20G×1-1/8"	0.9×28mm		Yellow
SSVS19/	5VS19	19G×1-1/8"	1.1×28mm		Cream
SSVS18/	5VS18	18G×1-1/8"	1.2×28mm		Pink

### Infusion therapy



### **Infusion set**







Delivering solutions in IV therapy for improved patient care.

#### Features

- ➤ Sharp spike is easy to puncture. ➤ The ISO standard spike can
- fit with any infusion bottles and bags.

  Soft and transparent
- ➤ Soft and transparent chamber is easy to press and adjust solution.
- ► Antibacterial air-inlet.



• Big chamber with air vent, filter, needleless Y site, Adelberg flow regulator, 1.8m tube, luer lock with cap.



 Big amber-colored drip chamber with air vent, ABS Adelberg flow regulator, amber-colored tube 1.5m, luer lock with needle.



 Small drip chamber with air vent, Y site, 1.5m tube, luer slip without needle.



 Winged spike, medium drip chamber with air vent, filter, ABS Adelberg flow regulator, Y injection port, 1.8m tube, luer lock without needle.



 Medium drip chamber without air vent, 60 drops, filter, Y injection port, 1.8m tube, luer slip without needle.

# Infusion therapy **Infusion set** Spike with drip chambers Injection ports We offer a wide variety of IV accessories to meet the needs of current practice: ➤ Small drip chamber 14.30\*50mm. Medium drip chamber 15.50\*54mm. Big drip chamber 17.20\*61mm. With air vent, without air vent. ► Transfer set. 15drops, 20drops, 60 drops/ml. With filter or without filter. ► Small flow regulator, Medium flow regulator, ABS Adelberg flow regulator. Accurate flow regulator rate of 5mL/hr to a maximum of 300mL/hr. Latex straight tube, injection bulb, Y injection port. Latex free straight tube, injection bulb, Y injection port. Needless Y injection port, Needle free injection cap. Luer slip, luer lock connector. ▶ Red cap, blue cap, white cap.▶ PVC tube length 1.3m to 2m or any required length.

► PVC tube DEHP FREE. Non-PVC tube.

► Customization mold construction for any kind of components.

### Infusion therapy







- ▶ Soft and transparent chamber for easy to read graduated scale.▶ Floating disk for easy fluid level indication and auto shut-off to
- resist the air bulb to the vein.

  Integral 15 micron air-guard filter to prevent air introduction.
- ► 60drops/ml. Gravity use.





Soft chamber 100ml, with air vent, latex tube, tube 1.5m, luer lock with needle.  Hard chamber 100ml, with air vent, latex bulb, luer slilp with needle  Soft chamber 110ml, with air vent, with latex shut-off valve, Y injection port, tube 1.8m, luer lock with needle  IVB 100-1  IVB 100-1	Description	Code
Soft chamber 110ml, with air vent, with latex shut-off valve, Y injection port, tube 1.8m, luer lock with needle IVB 110-1		
Soft chamber 150ml, with air vent, with latex free shut-off valve, needleless Y injection port, tube 2.2m, luer lock without needle. IVB 150-2	Soft chamber 150ml, with air vent, with latex shut-off valve, Diameter 5.6mm tube to comply with the infusion spike.	IVB 150-1

### Infusion therapy



### IV cannula







Delivering solutions in IV therapy for a short term period.





- ► Color-coded casing cap allows for easier identification of catheter size.
- ► Translucent catheter hub and flashback chamber allows for easy detection of blood.
- ► Flash back at vein insertion.
- Teflon radio-opaque catheter.
   Precision finished pet or pu catheter assures stable flow and eliminates catheter tip kink.
- ► The sharp needle features a universal bevel geometry that allows a wider choice of insertion
- ► Can be connected to syringe by removing filter cap to expose luer taper end.
- ▶ Use of hydrophobic membrane filter eliminates blood leakage.
- Close and smooth contact between catheter tip and inner needle enables safe and smooth.
- ▶ Vein puncture.

	Code		Gaug	e	Size and length	Color
			V			3
Pencil type	Wing type	Injection type				
IVCP14	IVCW14	IVCI14	14G		2"	Orange
IVCP16	IVCW16	IVCI16	16G		2", 21/2"	Grey
IVCP18	IVCW18	IVCI18	18G		11/4", 11/2", 2"	Green
IVCP20	IVCW20	IVCI20	20G		11/4", 2"	Pink
IVCP22	IVCW22	IVCI22	22G		1"	Blue
IVCP24	IVCW24	IVCI24	24G		3/4"	Yellow

### Infusion therapy



### Safety IV cannula







- ▶ Safety design can reduce needle stick Injuries.
  ▶ Result in less blood exposure to staff than catheters featuring spring-retraction of the needle.
- ▶ Promote first Stick Success.
- The safety design ensures that the safety mechanism is always activated and prevents re-insertion of the stylet.
- ▶ No extra steps are required to ensure healthcare worker and patient safety.
- ▶ With this easy-to-use safety IV catheter, you can be confident of successful needle and catheter placement every time.
- PVC-Free, DEHP-Free and Latex-Free.
   Close system and open system available.







Size and Length	Code
	TYPE 1 TYPE 2
18G×30mm, 18G×48mm	IVCS1-18 / IVCS2-18
20G×30mm, 20G×48mm	IVCS1-20 / IVCS2-20
22G×19mm, 22G×25mm	IVCS1-22 / IVCS2-22
24G×19mm	IVCS1-24 / IVCS2-24
26G×19mm	IVCS1-26 / IVCS2-26

### Infusion therapy



### IV cannula dressing







Fixing IV catheter, PICC and CVC catheter. Apply for friction skin or wet skin, such as skin care after tracheal incision. Apply for slight bruise, cut, suturing wound. etc. as a secondary



#### Features

- ► High permeability, avoid causing maceration to the skin.
- ▶ Prevent bacteria invasion, diminish the risk of infection.
- ▶ Waterproof, patient can take shower or bath.
- ► Transparent, wound situation could be observed anytime.
- ► Soft and comfortable, suit for different body parts, hypoallergenic.
- ► Strong viscosity, could keep on skin for 7days.
- Safe and convenient, the design of dressing make it easy operation
   Various sizes and styles, including paper frame IV dressing, new Utility style, hydrocolloid IV dressing, high-permeable grid-coating transparent dressing, PE foam frame IV dressing, non-woven IV dressing, and could fit to different clinical needs.
- ▶ Easy to use, could prevent wrinkling of dressings (absorbent pad could add to the dressing based on customer's requirements.
- ► Super absorbency, can prevent maceration.

#### Sterile IV start kit



#### Description

- ▶ 2ea Gauze sponges, 2x2"8ply.
- ▶ 1ea Label with four lines (optional).
- ► 1ea Tourniquet 3/4x18 synthetic.
- ▶ 1ea IV dressing film 6x7cm.
- ▶ 1ea PE tape, 3/4x18".
- ► 1ea Prep. appl 2% chg, 1,5ml (optional).
- Customization types are possible.

### Infusion therapy



#### **Extension tube**







- Our extension sets are available in various lengths and dimensions to accommodate
- ▶ A range of lengths that provide flexibility in a variety of clinical
- ► Available with various needleless connector options.
- A range of various flow clamper, components, connector options.
- ► To prevent dislodgment of IV catheter and offers ease of use during connection.
- ► Components are available for DEHP free or latex free.









Description	Code
Needleless type, 1way, female luer lock, male rotating luer lock, slide clamp,OD 4mm	ET01
Needleless type, 2way, female luer lock, male rotating luer lock, slide clamp,OD 4mm	ET02
Needleless type, 3way, female luer lock, male rotating luer lock, slide clamp,OD 4mm	ET03
Needleless type, 4way, female luer lock, male rotating luer lock, slide clamp,OD 4mm	ET04
Male and female luer lock	ET05
T connector, male and female luer lock, slide clamp	ET06
3way stopcock, male rotating luer lock	ET07
Needleless type, 1way, female luer lock, male rotating luer lock, slide clamp,OD 2.35mm	ET08
Needleless type, 2way, female luer lock, male rotating luer lock, slide clamp,OD 2.35mm	ET09
Needleless type, 3way, female luer lock, male rotating luer lock, slide clamp,OD 2.35mm	ET10



### Three way stopcock







Maximize fluid flow control and fluid and drug delivery efficiency.

#### **Features**

- ► We can offer a full range of stopcocks to assist you in controlling fluid flow through IV delivery
- ▶ Single packaged, as various options to meet all
- of your IV therapy needs.

  All options are liquid resistant with a transparent housing for clear view of the fluid path.









### Syringe and needle



### Disposable syringe with needle








Injection or harvesting of biological material, depending on choice of needle.



- ► Luer slip or luer lock nozzle, central or eccentric nozzle.
- ► High transparency barrel.
- ► Smooth consistent plunger action.
- Clear, legible and easy to read graduated scale.
- Minimal dead space.
- Triple contact rubber gasket and plunger-block to prevent accidental leakage during aspiration.
   Well rounded finger grips for safe handling.







### Disposable syringe with safety needle







#### Features

- ➤ Compatible with all standard hypodermic needles.
  ➤ Convenient to use, one hand activation, minimal
- training required.

  Economic, reduce the cost of safety injection.

  Triple safety mechanism to guarantee 100% safely.

- confirmation of safety activation.
- ► Click device guarantee no obstruction. during injection.
- ► Both luer slip and luer lock nozzle are available.





## Syringe and needle



Disposable syringe with safety needle





DSS1-2 DSS2-2 DSS1-5 DSS2-5 DSS1-10 DSS2-10

DSS1-100 DSS2-100



Syringe and needle



2-3ml 5ml 10ml

20ml 30ml 50-60m 100ml

SS0.5 SS1

SS3 SS5

### Syringe and needle



### Safety syringe







It is used in extraction of clinical liquid medicine and for injection, help eliminate accidental needle stick injury.

#### Features

- Available without needle or needles. Compatible with all standard hypodermic needles.
- ▶ Convenient to use, safety mechanism will be activated when the intended dose has been delivered.
- ▶ Pull back the plunger after injection, the locking ring secures the needle in the barrel of the syringe preventing the reuse of syringes.
- ▶ Special design can prevent the needle back into the barrel occasionally.
- No gap design guaranteed the complete injection without dead space.
- Exchangeable needles allow clinician to change needle after aspiration of medication.
- Visual confirmation of safety activation.
- ▶ Various range of syringes to meet the clinical requirements.



0.5ml 1ml

3ml 5ml 10ml

#### Using step

- ► Take off the needle protector.
- Aspiration the medication .
- ► The safety medication is automatically activated.
- After injection, pull back the plunger to the locking ring.
- ▶ Break the plunger.



### Syringe and needle



#### Safety needle







Syringe and needle

It is compatible with all syringes in clinical injection of liquid medicine, easy to use, help eliminate accidental needle stick injury.

- ▶ It is compatible with all standard syringes with luer lock or luer slip tip.
- Optimum visualization of injection site.
- After the injection, only one hand activation, the needle was covered completely which can immediately achieve the purpose of protection.

  Special safety lock on needle hub, triple safety mechanism to guarantee
- 100% safe protection.
- ► Effective function and simple design can reduce the cost of safety
- ► Easy and quick identification of the needle gauge by ISO colorful
- ► A wide range of sizes are available.



- ► Assemble the safety needle on the luer tip of sterilized
- Open the safety protector, which will open automatically under the elastic force.
- ► Aspirate the medication.
- ► After injection, apply a force on the safety protector by one hand until the safety protector covers the needle completely.



Code	Gauge	Size and length	Color
V .	V	V	V
SN1-29	0.33×13	29G×1/2"	Red
SN1-27	0.40×10-16	27G×3/8"-5/8"	Grey
SN1-26	0.45×10-16	26G×3/8"-5/8"	Brown
SN1-25	0.50×16-25	25G×5/8"-1"	Orange
SN1-24	0.55×25	24G×1"	Purple
SN1-23	0.60×25-38	23G×1"-1-1/2"	Blue
SN1-22	0.70×25-38	22G×1"-1-1/2"	Black
SN1-21	0.80×25-38	21G×1"-1-1/2"	Green
SN1-20	0.90×38	20G×1-1/2"	Yellow
SN1-19	1.10×38	19G×1-1/2"	Cream
SN1-18	1.20×38	18G×1-1/2"	Pink
SN1-16	1.60×38	16G×1-1/2"	White



Safety needle











- ► Latest innovative design with special safety shield.
- ▶ Protection to healthcare worker during transporting device and after use.
- Give much easy operation to nurses.
- Safely locked by a slight pushing up.Compatible with all kinds of syringes.



Code	OD and length(mm)	Gauge and length(inch)	Code	OD and length(mm)	Gauge and length(inch)
SN2-27	0.4×16	27G×5/8"	SN2-21	0.80×25	21G×1"
SN2-25	0.50×16-25	25G×5/8"-1"	SN2-20	0.90×25	20G×1"
SN2-23	0.60×25	23G×1"	SN2-18	1.20×25	18G×1"
SN2-22	0.70×25	22G×1"			

### Syringe and needle



### Disposable hypodermic needle







Intradermal, subcutaneous, intramuscular and intra-articular injection, harvesting of biological material depending on chosen gauge.

- ► High quality 304 stainless steel.
- ► Triple sharpened tip with anti-coring treatment maximize patient comfort.
- ▶ Both short bevel and long bevel are optional.
- ► Silicon coated to ease skin penetration and minimise pain.
- Large ID, high flow, both thin wall and regular wall are available.
- ► Semi-transparent needle-hub.
- International color coded hub by size for clear recognize.
- Customer special needle size as per required.



	Size	Color	Y				Leng	th					IC		Bevel	Angle	Code
	1	V	W									7	<u>/</u>	V	V	V	
OD (mm)	Gauge	Color&Code	4 1/6"	6 ′ 1/4″	8 1/3"	13 1/2"	16 5/8"	20 3/4"	25 1"	32 1 1/4"			RW	TW	Long Bevel	Short Bevel	
0.25	31G	Blue	√	√	$\checkmark$	$\checkmark$	$\checkmark$						$\checkmark$		$\checkmark$	$\checkmark$	HN31
0.30	30G	Light yellow		$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$						$\checkmark$		$\checkmark$	$\checkmark$	HN30
0.33	29G	Red	√	√	√	$\checkmark$	$\checkmark$						$\checkmark$		$\checkmark$	$\checkmark$	HN29
0.36	28G	Blue-green	√	√	√	√	√						$\checkmark$		$\checkmark$	$\checkmark$	HN28
0.40	27G	Grey	√	√	√	$\checkmark$	√						$\checkmark$		$\checkmark$	$\checkmark$	HN27
0.45	26G	Browm				$\checkmark$	√						$\checkmark$		$\checkmark$	$\checkmark$	HN26
0.50	25G	Orange					√	$\checkmark$	√		$\checkmark$		$\checkmark$	√	$\checkmark$	$\checkmark$	HN25
0.55	24G	Purple						$\checkmark$					$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	HN24
0.60	23G	Blue							√	$\checkmark$	$\checkmark$			$\checkmark$	$\checkmark$		HN23
0.70	22G	Black							√	√	$\checkmark$			√	$\checkmark$		HN22
0.80	21G	Green							√	√	$\checkmark$	√		√	$\checkmark$		HN21
0.90	20G	Yellow									$\checkmark$	√		$\checkmark$	$\checkmark$		HN20
1.10	19G	Cream									√	√		√	$\checkmark$	$\checkmark$	HN19
1.20	18G	Pink									$\checkmark$	√		√	√	√	HN18
1.40	17G	Red-purple									√			$\checkmark$	$\checkmark$	$\checkmark$	HN17
1.60	16G	White									$\checkmark$	√		$\checkmark$	$\checkmark$	√	HN16
1.80	15G	Blue-gray									$\checkmark$	√		√		$\checkmark$	HN15
2.10	14G	Pale-green									$\checkmark$	$\checkmark$		$\checkmark$		$\checkmark$	HN14
Remai	rk: We can m	ake special size acc	cordin	ng to c	lients	' reque	est.										



# Auto-Disable syringe for fixed-dose immunization







It is mostly used for WHO or UN aid project, the plunger of syringe was locked to prevent the reuse, help eliminate touch contamination.

#### Features

- ▶ In all cases, once the auto-disable feature has been activated, It was not possible to re-use the syringe and the needle under the normal conditions of use.
- ▶ It was not be possible to override the auto-disable feature when tested in accordance with the test method of ISO7886-3 Annex C.
- ▶ If an integrated needle is used, it shall become an integral part of the syringe and can not be detached.
- ▶ The auto-disable feature is automatically and remains effective from the point when 50% of the intended fixed dose has been delivered.
- ▶ Invention patented.



Code	Size	Needle size
ADS1-0.05	0.05ml	23G-27G
ADS1-0.1	0.1ml	23G-27G
ADS1-0.2	0.2ml	23G-27G
ADS1-0.3	0.3ml	23G-27G
ADS1-0.4	0.4ml	23G-27G
ADS1-0.5	0.5ml	23G-27G
ADS1-1	1ml	23G-27G



### Syringe and needle



# Auto-Disable syringe for fixed-dose immunization







➤ Syringe with non-integrated needle shall have a minimum needle unicon force applied as push or pull in the direction of the needle axis in accordance with ISO 7864:1993.

▶ The auto-disable feature is automatically activated on completion of the delivery of the intended fixed dose.

► Central luer slip or eccentric luer slip.



Type 2, With non-integrated needle

Code	Size	Code No.	Size
ADS2-2	2ml	ADS2-5	5ml
ADS2-3	3ml	ADS2-10	10ml

### Pre-filled flush syringe

Our flush syringes provide dependable, cost-effective alternatives to vial-based flushing systems.







### Dental needle









- ▶ High quality 304 stainless steel.▶ Short and long needle bevel are
- ► Color coded hub by size for clear recognition.
- ▶ Different sizes available.
- ▶ Designed special for dental use.





### **Irrigation needle**

The irrigation needle are intended for irrigation purposes in dentistry and ophthalmology, or used for withdrawal of drugs/solutions out of vials.



### Syringe and needle



### Suture removal Set










Qtty.	Description					
1ea √	Metal Littauer scissors					
1ea	Plastic suture forceps					
1ea	Gauze pad					
1ea	Alcohol prep pad					
1ea	lodine prep pad					
Customization types are possible.						



### Veterinary needle

Used for veterinary intradermal, subcutaneous, intramuscular and intra-articular injection, harvesting of biological material depending on chosen

- ▶ Luer fitting or record fitting are available.
- ▶ By copper base and strong enough.
- Full range for chosen.







#### Hard pack syringe







We offer syringes for the pet care market.

#### Features

- ► Hard pack design, convenient and safe cartridge pack, easily disposed of and transportable.
- ▶ Heat sealed cap for assurance of sterility.

Code	Size
HS01	: 1ml
HS02	2ml
HS05	5ml
HS10	10ml



### Hazardous drug safety



#### Drug preparation







We offer various drug preparation needles for safety drug admixture, minimizes the risk of needle stick injury while drawing up medication, medication re-constitution and accessing IV bags.

#### Aspiration needle with safety shield

#### Features

- ▶ USA Patent pending, PCT with a priority date already initiated.
- ► FDA 510K approved.
- ► The shield that prevents accidental needle stick injury, 100% to protect the doctors and nurses from the danger.
- Intended for aspiration of medications from all kinds of vials safely.
- ► Injection of medications into IV rubber ports safely.
- ➤ Transfer of biological specimens, blood into vacuum specimen tubes safely.
- ► Compatible with standard luer slip and luer lock syringes.



FDA 510k approved

#### Blunt needle with/without filter

Blunt needle is a sterile, single use device used for safety drug admixture.

- ▶ 45° bevel angle optimized for fluid transfer with e.g. stoppered vials, glass and plastic ampoules.
- ▶ 45° bevel angle thin walled cannula for faster draw-up (up to 20% higher draw-up speed in comparison to comparable products in the market).
- ▶ Made of stainless steel with special smoothing bevel.
- Blunt smoothing bevel design minimizes the risk of needle stick injury while drawing up medication, medication reconstitution and accessing IV bags.
- ► Red colored hub and protective cap for clear differentiation.
- ▶ Blunt needle with filter with a 5micro filter membrane to prevent glass particles being drawn into the syringe when preparing medication from ampules.



### Hazardous drug safety



### Drug preparation









#### **Curved needle**

#### Features

- ▶ The curved needle will reduce the risk of needle stick injuries.
- ▶ To avoid cutting the rubber pieces during the penetration.



#### Side hole needle

#### Features

► The side hole needle could avoid cutting the rubber pieces during the penetration.

Spec.	Size	type	code	
. Blunt needle	: 14G	: without filter :	BLN14	-
Blunt needle	16G	without filter	BLN16	
Blunt needle	18G	with filter / without filter	BLNF18/BLN18	
Blunt needle	19G	with filter / without filter	BLNF19/BLN19	
Blunt needle	20G	without filter	BLN20	
Blunt needle	21G	without filter	BLN21	
Curved needle	14G		CVN14	
Curved needle	15G		CVN15	
Curved needle	16G		CVN16	
Curved needle	17G	i i	CVN17	
Curved needle	18G		CVN18	
Side hole needle	16G		SHN16	
Side hole needle	18G		SHN18	
Aspiration needle	18G		ASN18	

### Hazardous drug safety



### Drug preparation















- ▶ Needle-free system that utilizes split septum technology similar to that of a conventional needle system.
- ► A blunt plastic cannula replaces the standard needle, to provide needle-free safety while requiring little or no change in existing protocol or IV procedures.
- ▶ Designed to eliminate accidental needle stick.
- ▶ Provides a swabbable injection site that is a closed system, accessed by inserting a cannula.
- Designed to ensure a safe, secure connection with Clip Lock cannula.
- ► Compatible with IV fluids, including blood and lipids.
- Visually inspect with a clear outer housing that allowsvisible flush confirmation, and an open fluid path that hasno fluid entrapment areas.
- ► Clean its smooth access surface facilitates effective cleaning and disinfection.
- ► Connect minimal connection force is needed, allowing quick, effortless access and easy syringe connection.
- ► We offer a variety of needleless connector devices to meet the desired need of individual facilities.



Code	Description
Safesite	Compatible with IV fluids, including blood and lipids.
SafesiteY	Compatible with infusion set and IV burettes
SafesiteP	Innovative positive displacement feature is designed to help reduce catheter occlusions while preventing healthcare workers from accidental needlestick injuries. Compatible with IV catheter.
Safespike	Compatible with all infusion bottle or medication vial, the ergonomic grip design is ideal for clinician handling and patient comfort, offering an easy-grip barrel designed to minimize slips and touch contamination, help make reconstituting, admixing and IV drug delivery secure and cost-effective.
SafespikeF	With filter, help make reconstituting, admixing and IV drug delivery secure and cost-effective. Help ensure fluid integrity by removing contaminated particles with a broad range of filtration devices.  1, Min. dia. of filter: 0.45µm  2, Min. dia. of filter: 5µm

### Hazardous drug safety



#### **Sharps disposal solutions**









# Safety box for the disposal of used sharps

It could contain used contaminated sharps (including syringe, needles, phlebotomy devices, suture needles, IV insert needles, butterflies, lancets, scalpels) safely and efficiently, until safe disposal or recycling:

-At the point of use

-During temporary storage

-During handling and transport to the point of treatment and final disposal.

#### Features

- ► Exactly meet the requirements of WHO/PQS/E10/SB01.1.
- ▶ Resistance average force of penetration more than 15N from each position.
- ▶ Resistance damage during 100 drops from height of 0.8m.
- Resistance temperatures of up to 170°C periods up to 30minutes.
- ▶ Resistance to 48 hours at 43 °C and 90% relative humidity in 5mm of water.
- ▶ Could be supplied flat-packed to minimize shipping and storage volume.
- Clearly marked with international bio-hazard warning.
- Could be carried in one hand without spillage of contents and without risk of needle stick injury.
- ► The Bio-degradable cardboard recycled material is used.
- Non-toxic inks, glues and dyes.
- ► Clear pictorial instructions without writing.
- Chinese invention patented.

Code	Size	Sharp aperture dimension
SBS1.25	1.25L	: Dia. 38mm
SBS2.5	2.5L	Dia. 38mm
SBS5	5L	Dia. 38mm
SBS10	10L	Dia. 38mm
SBS15	15L	Dia. 38mm
SBS20	20L	Dia. 38mm

#### • Colors: Unbleached sulphate board, or non-chlorine bleached white, or yellow

## Hazardous drug safety



#### Sharps disposal solutions







#### **Sharp container**

The waste containers used for disposal of hospital waste, it is characterized by an elevated level of safety and reliability. The material used for production guarantee the solidity, the hermetic closure, and the high resistance to shock and puncture. The containers are equipped with disposable lids.







-		<u></u>	
Code	Volume (L)	Diameter (cm)	Height (cm)
SC-1R	1.0	11.5	12.5
SC-2R	2.0	13.5	14.5
SC-3R	3.0	16.5	17
SC-4R	4.0	17.5	18.5
SC-5R	5.0	19.5	20
SC-6.5R	6.5	21.0	21
SC-8R	8.0	23.0	23
SC-15R	15.0	27.0	27
SC-3.5S	3.5	20×16×13	
SC-5S	5.0	25×20×15	
SC-8S	8.0	30×22×16	

### Hazardous drug safety



#### **Sharps disposal solutions**







#### **Ultimate sharps safety stations (needle counter)**

It is used for needle or sutures counting in the operation room, dangerous sticks from uncounted needles. Can be safely avoided from the needle stick injury.



#### Features

- ▶ We offer a wide selection of counting and disposal systems that utilize foam block, foam strip, adhesive or magnetic retention.
- ► Each system is housed in a flexible, impactresistant plastic box that provide protection while minimizing unnecessary handling.
- ▶ We have USA PETENT for needle counter with more safety and function, exclusive authorized manufacturer in China.
- ▶ PCT with a priority date was already initiated.
- ▶ Detachable halves allow the same needle counter to be used in two separate areas.
- ▶ 2 universal needle safety holding device enable to hold either disposable syringe or scalpel.
- ► Surgical blade safety removal device with a protection wall to protect the iumping of blades.
- ► A needle sheathing safety device to help user's hand free and clear of the needle.
- ► The safety lock device guarantees the needle or blade will not slip from the gap of needle counter.
- ▶ The combination of latch and hinge provides a surprisingly secure, simple and inexpensive integrally safety box.



Code	Counts	Spec.
		V
NC10	10	
NC15	15	Magnet/ foam block/
NC20	20	foam strips/adhesive
: NC30	30	: magnet strip/ clear lid
NC40	40	are optional
NC60	60	
: NC80	80	
UNC20-1	20	Ultimate with safety
: UNC20-2	20	: needle, blade removal
UNC40-1	40	device, and safety
UNC40-2	40	lock device

#### **Endoscope transport pad**

It is used to pack the disposal medical drape and kits after surgery, with the drawstring design which is easy for disposal.

Dimension: 93\*100cm with core size: 41\*100cm, Paper + PE material 90g.

Any customization sizes are available.



### Hazardous drug safety



#### Safety medical products for **NEEDLE STICK PREVENTION**







Although the Needle Stick Safety and Prevention Act was signed into legislation in the United States in 2000, there still remains much to be done to protect health care workers against sharps injuries that potentially expose them to deadly blood borne pathogens. According to recent research, (2011 Healthcare Worker Risks from Short Peripheral Catheter Insertion), percutaneous injuries to healthcare workers from a sharps device are four times more likely to occur than mucocutaneous exposure and pose a greater threat to disease transmission. Also, according to the 2008 Study of Nurses Views on Workplace Safety and Needle Stick Injuries, nearly two-thirds of US nurses say needle stick injuries and blood borne infections remain major concerned.

In order to protect health care workers against sharps injuries, KB Medical Group is committed to address these issues by providing the most effective safety-engineered therapy devices for optimum safety in the hospitals.

#### Our safety family products includes: (For detail product information, go to the page mentioned below):

- ► Safety scalp vein set (Page 1).
- ► Safety IV cannula (Page 6).
- ► Safety syringe (Page 11-13).
- ► Safety needle (Page 14-15). ► Auto-Disable syringe(Page 17-18).
- ► Safety drug preparation needle (Page 22-23).
- Needle free system (Page 24).
- ► Safety box for the disposal of used sharps (Page 25).
- ► Sharp container (Page 26,52).
- ▶ Ultimate sharps safety stations (Page 27).
- ► Safety blood collection needle (Page 34).
- Safety blood lancet (Page 35).
- ► Safety fistula needle (Page 47).
- ► Safety insulin syringe (Page 49-50).
- ► Safety insulin pen needle (Page 51).





# KB

IS080369-3









Low dose

With the further research on the structure and function of gastrointestinal tract, it has been gradually recognized that gastrointestinal tract is not only an organ of digestion and absorption, but also an important immune organ.

Enteral feeding refers to the delivery of a nutritionally complete feed, containing protein, carbohydrate, fat, water, minerals and vitamins, can directly into the stomach, duodenum or jejunum.

Enteral nutrition is the nutritional support provided by the gastrointestinal tract to provide nutrients needed for metabolism and other nutrients. It depends on the length of time, mental state and gastrointestinal function. There are three ways for enteral access:

- 1, Nasogastic feeding.
- 2, Gastrostomy.
- 3, Jejunostomy.

We offer the various components for Enteral nutrition therapy.

#### **ENFit®** syringe

#### Features

- ► The Enteral ENFit® syringe tip is compliant with the ISO80369-3 to help reduce the risk of misconnections.
  - Offer a full range and are the result of research aimed at creating a product with certified quality and safety features.
    - ▶ DEHP Free.
    - ► Latex free plungers, made of synthetic rubber which does not contain natural latex, in order to prevent any allergies.
      - ► Low dose syringes have a male design feature within the ENFit® female connector to prevent volume in the dead space.
        - ► Tip cap is available. (ENC6, Page30).

	Spec.	Size	Code
	Enfit and low dose	1ml	ES1
	Enfit and low dose	3ml	ES3
-	Enfit and low dose	5ml	ES5
	Enfit	10ml	ES10
	Enfit	20ml	ES20
	Enfit	30ml	ES30
	Enfit	60ml	ES60
	Bulb	60ml	BS60
	Thumb control ring	60ml	TS60
1	Catheter tip	60ml	CT60

## Nutrition therapy



### **ISO80369-3**

#### CE





#### The ENFit® connector

#### Features

- ► The ENFit® connector is compliant with the ISO80369-3 to help reduce the risk of misconnections.
- ► The ENFit® connector make enteral feeding and medication administration easy, efficient and most importantly, comfortable for patients.
- Have a wide range of connector for enteral feeding which along with the syringes, complete a range of products which enable solutes and medication to be prepared and administered in complete safety.
- Customized different connectors according to requirements.











30

on therapy









### Nutrition therapy



#### **ISO80369-3**

#### **ENFit®** Feeding set

#### Features

- ► ENFit® feeding sets are compliant with the ISO 80369-3 to help reduce the risk of misconnections.
- ▶ Patient label provides convenient patient documentation.
- Latex-free and DEHP-free.
- ► EVA material available.

#### Advantage of EVA material

- ► Ensure absolute Non-DEHP, Non-TOTM, Non-DOA or any other
- Ensure absolute non risk of toxicity.
- ► Ensure the maximum safety for enteral feeding technology.



- ► Feeding bag (1000/1200/1500ml).
- ▶ PVC tube 1.73m /2.28m / customized.
- ▶ Regulator with wheel.
- ▶ Drip chamber.
- Y connector with cap. ► ENFit® lock female

11cm / 13cm / customized.

- connector with cap. ▶ Patient label.
- ► Silicon insert( optional )



- ► PVC tube 1.73m /2.28m /

- ➤ Silicon insert( optional ) 11cm / 13cm / customized.
- ▶ Feedcap adaptor.
- customized. ▶ Regulator with wheel.
- ▶ Drip chamber.
- Y connector with cap. ► ENFit® lock female.
- connector with cap.
- ▶ Patient label.
- ► Spike with cap.

► ENPLUS spike with cap.

▶ Regulator with wheel. ▶ Drip chamber. Y connector with cap.

▶ Patient label.

customized.

▶ PVC tube 1.73m /2.28m /customized.

► ENFit® lock female connector with cap.

► Silicon insert( optional ) 11cm / 13cm /

- PVC tube 1.5m/customized.
- ▶ Drip chamber.
- Regulator with wheel.
- ▶ Feeding connector with cap.



### Nutrition therapy



# IS080369-3

#### **Oral syringe**

#### Features

- ▶ We can offer both clear and amber syringe.
- ► Special tip design for wrong route administration prevention.
- Clear dose scale for accurate volume confirmation.
- ► Amber barrel design for protect light-sensitive
- ► Tip cap is available.



Specification	Size	Code
		Type 1 Type 2
Oral tip	1ml	ORS1-01/ORS2-01
Oral tip	3ml	: ORS1-03/ORS2-03
Oral tip	5ml	ORS1-05/ORS2-05
Oral tip	10ml	ORS1-10/ORS2-10
Oral tip	20ml	ORS1-20/ORS2-20
Oral tip	30ml	ORS1-30/ORS2-30
Oral tip	50ml	ORS1-50/ORS2-50



#### **Irrigation syringe**

#### Features

- ► Irrigation syringes are designed for simplicity and efficiency of care.
- ► Irrigation syringes are available sterile or non-sterile to suit your clinical need.
- ▶ Includes syringe tip protector.

For detail product information go to page 46.






## Specimen collection

## Blood collection needle











▶ Double wings make puncture more safe.

▶ Different color wings means different sizes.

► Soft and transparent tube, can observe the vein blood flow clearly.



- ▶ Pen type blood collection, latex free.
- Multi-sample needles permit several samples to be taken with a single puncture.
- ► Sharp & smooth edges make penetration painless, easy connection to rubber stoppers.



▶ Blood collection needle with hypodermic needle holder, latex free.

## Specimen collection



#### Safety blood collection needle








34

Infusion therapy

It is used for aspiration of blood, the safety mechanism help to eliminate accidental needle stick injuries.

- Type 1 Scalp vein set type.Type 2 Safety needle type.Type 3 Pen type.





Size		Code						
V	V	V 1	Nomal b	Nomal blood collection needle		Safety b	lood collection	n needle
OD(mm)	Gauge	Color	Type1	Type2	Type3	Type1	Type2	Type3
0.40	27G	Grey	BCN27-1	BCN27-2	BCN27-3	SBCN27-1	SBCN27-2	SBCN27-3
0.45	26G	Brown	BCN26-1	BCN26-2	BCN26-3	SBCN26-1	SBCN26-2	SBCN26-3
0.50	25G	Orange	BCN25-1	BCN25-2	BCN25-3	SBCN25-1	SBCN25-2	SBCN25-3
0.55	24G	Purple	BCN24-1	BCN24-2	BCN24-3	SBCN24-1	SBCN24-2	SBCN24-3
0.60	23G	Blue	BCN23-1	BCN23-2	BCN23-3	SBCN23-1	SBCN23-2	SBCN23-3
0.70	22G	Black	BCN22-1	BCN22-2	BCN22-3	SBCN22-1	SBCN22-2	SBCN22-3
0.80	21G	Green	BCN21-1	BCN21-2	BCN21-3	SBCN21-1	SBCN21-2	SBCN21-3
0.90	20G	Yellow	BCN20-1	BCN20-2	BCN20-3	SBCN20-1	SBCN20-2	SBCN20-3
1.10	19G	Cream	BCN19-1	BCN19-2	BCN19-3	SBCN19-1	SBCN19-2	SBCN19-3
1.20	18G	Pink	BCN18-1	BCN18-2	BCN18-3	SBCN18-1	SBCN18-2	SBCN18-3

Remark: normal size in above table, we can also make special size as clients' request.

### Specimen collection



Safety blood lancet









It is used in blood collection sampling, and help eliminate accidental needlestick injuries.

#### **Product Features**

- ▶ Needle is fully shielded before and after use.
- Cannot be reused.
- ▶ Sterilized by gamma ray.
- ▶ Pre-loaded operating design ensures a good blood sample.
- Ultra sharp tri-bevel needle, high speed penetration.
- ► Multiple depth selection with color coding.



	Code	Code	Size	Length
TYF	PE 1	TYPE 2	V	V
SBLC	)1-28	SBL02-28	28G	1.8mm
SBLC	)1-26	SBL02-26	26G	1.8mm
SBLC	)1-21	SBL02-21	21G	1.8/2.2/2.8mm
SBLC	)1-18	SBL02-18	18G	1.8mm

New item, less pain
Diabetes clinics/hospital wards for blood glucose monitoring
Haemoglobin, cholesterol blood typing, coagulation, blood gas and other diagnostic tests, For paediatrics and normal adult skin
Special baby care and maternity units for heel puncture-neonates and preemie

Application

### Specimen collection



### **Blood collection tube**







36

Infusion therapy

We provide a full offering of blood collection tubes in both glass and plastic.

Blood collection tubes include silicone coated tubes that indicate the interior surface of the blood collection tube has been coated with a water soluble silicone material.



Color	Additive	Application	Volume	Size	Color	Additive	Application	Volume	Size
Red	Plain	Serum blochemistry	3ml	13×75	Green		Plasma blochemistry	2ml	13×75
		Electrolyte Tumor mark		13×100		Blood gas analysis Blood stream changing		13×100	
		Medicament inspecting						3ml	
		Thyroid function						4ml	
		AIDS inspecting						5ml	
			5ml		Purple	EDTA-K2	Blood routine	1ml	13×75
			7ml			EDTA-K3	Examination(Red blood		13×100
Red	Coagulant	Serum blochemistry	3ml	13×75		EDTA-Na2	cell.White blood cell.		
		Electrolyte		13×100			Blood platele, Classify of		
		Tumor mark					blood cell)		
		Medicament inspecting		:				2ml	
		Thyroid function						3ml	
		AIDS inspecting		:				4ml	
			5ml			- 00 ( 5 1)	=	5ml	
\ / II			7ml		Black	3.8% Sodium	, ,	2ml	13×75
Yellow	Separating gel	Serum blochemistry	2ml	13×75		Cltrate	sedimentation rate(ESR)		13×100
	1.0	Electrolyte		13×100	DI	- 00/ 5	DT 4 DTT TT 5'	2.4ml	
		Tumor mark			Blue	3.8% Sodium	,	1.8ml	13×75
	activator	Medicament inspecting,					Platelet coglutionation		13×100
		Thyroid function				solution	Testing	27	
		AIDS inspecting	71					2.7ml	
			4ml		DI .	7.00/.6.1	District and Constitution	3.6ml	4775
\\/\bito	6 11	Dland sugar tast	6ml	: 17975	Blue		Platelet coglutionation	1.8ml	13×75
White	: Sodlum : floride/	Blood sugar test	3ml	13×75		Cltrate buffer	resting(PAGT)		13×100
	potassium	Tolerence tasting for blood		: 13×100		solution		27	
	oxalate	: sugar	:	:		:		2.7ml	:

## 



#### Spinal needle

Subdural, lower thorax and lumbar spinal puncture.

#### Features

- ▶ Full sizes of anesthesia needles.
- ▶ Identified as spinal needle bevel quincke tip, pencil point tip and epidural needle.
- ▶ 304 stainless steel needle tubes and stylets.
  International colour coding.
- ► With mounted introducer needle.
- Luer lock hub.





Quincke tip	Introducer	Length	Code
18G×3-1/2"	None	38-220mm	SNQ18
19G×3-1/2"	None	38-220mm	SNQ19
20G×3-1/2"	None	38-220mm	SNQ20
21G×3-1/2"	None	38-220mm	SNQ21
22G×3-1/2"	None	38-220mm	SNQ22
23G×3-1/2"	19G×1-1/2"	38-220mm	SNQ23
24G×3-1/2"	20G×1-1/2"	38-220mm	SNQ24
25G×3-1/2"	20G×1-1/2"	38-220mm	SNQ25
26G×3-1/2"	20G×1-1/2"	38-220mm	SNQ26
27G×3-1/2"	22G×1-1/2"	38-220mm	SNQ27

- ► Trauma-free OUINCKE short tip to ensure immediate injection point closure and reduced CSF loss.
- ▶ Differentiated angle coupling of needle and obturator tip to avoid production and movement of tissue fragments, improving operator



Pencil point	Introducer	Length	Code
22g×3-1/2"	None	90-220mm	SNP22
23G×3-1/2"	19G×1-1/2"	90-220mm	SNP23
24G×3-1/2"	20G×1-1/2"	90-220mm	SNP24
25G×3-1/2"	20G×1-1/2"	90-220mm	SNP25
26G×3-1/2"	20G×1-1/2"	90-220mm	SNP26
27G×3-1/2"	22G×1-1/2"	90-220mm	SNP27

- ▶ With no sharp edges separates fibers without slicing.
- ▶ The lateral hole shape and its position relative to the tip favour rapid CSF reflux and effective analgesic administration.
- ▶ Reduce the risk of saddle block.

## Anesthesia & pain management (€

#### Anesthesia kit







Single-use anesthesia set is clinically used for epidural or spinal anesthesia during puncture and drug injection. This product can only be for single-use, in order to avoid cross-infection. It contains the unique components needed to perform anesthesia operation.



#### Components

- ▶ Epidural needle, Size: 15G,16G,17G,18G. Length: 80mm, 90mm.

- ► Epidural filter: 0.22µm with hydrophobic membrane luer lock.
- LOR syringe: 5ml, 7ml, 10ml.
- Catheter adaptor: TBA, secure the catheter completely.

#### Optional

- ▶ Epidural kit, spinal kit and combined kit available.
- Catheter with /without metal wire available.
- ► Precision filter 0.2-5µm available.
- ▶ Double layer UV filter 0.2-5µm are available.
- Any customization components.



# Anesthesia & pain management (Fig. 4)




### High pressure manifold









### NRFit™ syringe

It is used for regional anesthesia, the new patient safety initiative neuraxial products with NRFit™ connectors according to American Society of Anesthesiologists (ASA) meeting, Oct. 22-26 in Chicago 2016. All healthcare facilities are expected to begin converting to meet the ISO standard requiring NRFit connectors for neuraxial use. Healthcare facilities has already taken legislative action by mandating that convert to these new standardized epidural products by Jan. 1, 2017.



- ► NRFit<sup>™</sup> connectors design--ISO80369-6 compliant.
- ► Neuraxial applicant connector.
- ► Small bore connectors that are intended to prevent misconnections to improve patient safety.

#### NRS3 NRFit™ connector NRFit<sup>™</sup> connector NRS10 NRFit™ connector NRFit™ connector

### Transfusion therapy



### **Gravity blood administration sets**







Filtration for various blood IV administration needs.

#### Features

- ▶ Sharp spike is easy to puncture.
- ▶ Plastic spike/metallic spike available.
- ▶ Vented with hydrophobic bacterial filter.
- ► Soft and transparent chamber.
- ► Extra soft & kink resistant tubes.
- ► Acute roller clamp.
- ▶ With 200micron filtration.
- Latex free.
- ► Straight and Y type are available.
- ▶ DEHP free tube is available.



Big chamber with air inlet



Big chamber without air inlet



### Transfusion therapy



#### **Blood bag**

It is used to collect and preserve whole human blood and blood components in hospital or blood bank.









#### Features

- ▶ Our blood bag could not only guarantee high quality of the blood components but also to provide maximum safety during blood donation and sample collection.
  - ► The safety needle shield enhances safety against the accidental needle stick injury that can transmit HIV, Hepatitis B and C and other blood borne diseases.
    - ► Closed system & reusable blood sampling with vacuum tube.
    - ► Label material with precutting line and could not be tear integrated.



Single blood bag is used for collection of whole human blood. It consists of one bag containing anticoagulant solution CPDA.

#### BB02/Double

Double bag system consists of one primary bag containing anticoagulant solution CPD/CPDA and one empty satellite bag, for separation of half plasma.

Triple bag system consists of one primary bag

anticoagulant solution CPD/CPDA, one satellite bag containing SAGM preservation solution, and one

satellite bag. This additive can maintain red blood cell viability, with over 80% of survival rate after 5 weeks of storage.

#### BB04/Quadruple

Quadruple bag system consists of one primary bag containing anticoagulant solution CPD/CPDA, one satellite bag containing SAGM preservation solution, and two empty satellite bags. The additive SAGM can maintain red blood cell viability, with over 80% of survival rate after 5 weeks of storage.



### Transfusion therapy




#### Cord blood bag







It is designed to collect the umbilical cord blood (also known as placental blood) after the birth of the baby, within the shortest time and minimal risk of contamination. The clinical use of the cord blood is allergenic transplantation in patients suffering from severe blood diseases such as leukemia, lymphoma, etc.



#### Optional

- ▶ Top and top system.
- ▶ Top and bottom system.
- ► Empty transfer bag.
- ▶ PALL leukocyte filtration system.
- ▶ Blood sampling with vacuum tube.
- ► Customization bags.

With pull outlet

2000ml

UDBRP20

### Urinary care



### Urine drainage bag









- ► Anti-reflux drip chamber on the bag helps prevent urine reflux.
- ▶ With airvent and non-return valve ► Varies choice of components:
- universal hanger, bed sheet clamp, dog house, needless samples port.
- ► High grade PVC.
- Latex free.
- ► Tyvek pouch for O.R. pack or PE

UDBA20 ● Apple shape 2000ml



without outlet

2000ml

### Urinary care



### Leg bag



Rectangle or arc style for customer option.
It is secured with two comport straps in latex free, they have a twist drainage outlet or flip flow valve for easy emptying.

T-tap outlet or twist outlet for one hand easy emptying.

500ml

600ml

750ml

900ml

1000ml

With white fabric for comfort skin touch.

► Tube length 30-90cm for customer option.



### **Urine meter**

#### Features

Code

UM01

UM02

UM03

LB600

LB750

LB900

LB1000

▶ Designed for easy bag suspension with a reinforced double hanger and built -in anti-back

Spec.

200ml meter with 2000ml bag, 1 chamber. 500ml meter with 2000ml bag, 1 chamber.

Remark:200ml-500ml capacity meter is available

- ➤ Convenient needleless sample port is accessible with either needleless or blunt cannula, and offers added protection for patient and caregiver.
- Adjustable bed sheet clamp secures
- ► Include additional rope hanger.



### Urinary care



### Irrigation tray







Irrigation tray is designed to improve simplicity and efficiency of care.



#### Components

Qtty	Description	Qtty	Description
1ea	1000ml collection tray	2ea	Alcohol prep pad
1ea	500 ml graduated container	1ea	Moisture-proof under pad 45x45cm
1ea	60ml bulb or thumb control ring syringe	1pr	Vinyl gloves
1ea	Protector cap		, ,

### Urinary care



### Irrigation syringe







46

#### **Bulb syringe**

#### Features

- ▶ Flexible soft bulb promotes one handed administration.
- Designed to stand on the bulb end to help minimize contamination.
- ▶ Patient label provides convenient patient documentation.
- ► Sterile, Latex free, DEHP free, and BPA free.
- ► Includes syringe protector cap.

#### Thumb control ring syringe

#### Features

- ► Thumb control ring top allows for maximum control and one handed administration.
- Precision crafted O-ring plunger ensures a positive seal
   Re-sealable pole bag designed for safe and convenient
- ▶ Patient label provides convenient patient documentation.
- ▶ Sterile, Latex free, DEHP free, and BPA free.
- ▶ Includes syringe tip adapter and protector cap.

#### Flat top piston syringe

- ▶ Precision crafted O-ring plunger
- ensures a positive seal.

  Visible dosing marks designed for ease of use.
- ▶ Latex free, DEHP free, and BPA free.
   ▶ Includes syringe luer tip adapter and protector cap.

Size Description			Code
60ml	Thumb control ring type	V	URTS60
60ml	Bulb type		URBS60
60ml	Flat top piston		URFS60
12ml	: Curved tip	1	URCS12



## Renal therapy




### Fistula needle









We offer innovative, high performance products and services for use in dialysis. The core of the Renal Therapies spectrum provides a basis which targets improved clinical outcomes for the patient, while minimizing costs and risks for the healthcare providers.

#### Features

- Fixed and rotary wings are available.
  Japanese made sharp needle proves to be high quality.
  The large oval orifice on the needle back.
  Unique design and smooth finish.
  Super silicon coating of inner and outer needle wall to easy of penetration and minimize the pain.
- ► The spiral connector fits perfectly with dialysis line.

### Safety fistula needle

It is used for kidney disease patients to purify the blood, the safety mechanism help to eliminate accidental needle stick injuries.



	Code	Gauge	Size, length and color	Description
1	Non-safety/Safety			
-	AVFR15/AVFRS15	15G×1"	1.8×25mm Blue	Rotary wing with back eye
. <i>I</i>	AVFR16/AVFRS16	16G×1"	1.6×25mm Green	Rotary wing with back eye
-	AVFR17/AVFRS17	17G×1"	1.4×25mm Orange	Rotary wing with back eye
	AVFF15/AVFFS15	15G×1"	1.8×25mm Blue	Fixed wing with back eye
1	AVFF16/AVFFS16	16G×1"	1.6×25mm Green	Fixed wing with back eye
. /	AVFF17/AVFFS17	17G×1"	1.4×25mm Orange	Fixed wing with back eye

### Renal therapy



### **Blood line**





#### **Features**

- ► More than 500 types can fit with Fresenius, Cambro, B Braun, Nipro, Baxter hemodialysis machine.

  Luer type is tight with the adapter of fistula
- ▶ Different clamps fit for tubes of different
- Fixed connector with the wing is easy to join the adapter of hemodialysis machine.
- ► Customization model are possible upon the requirements.



Hemodialysis care kit



Qtty	Description	Qtty	Description
V	Before compartment	V	: After compartment
1ea	Surgical drape 45×45cm	1ea	Surgical drape 45×45cm
1pr	Examination glove	1pr	Examination gloves
5ea	Gauze ball plain type 4"×4"	2ea	Gauze ball plain type 4"×4"
1ea	Alcohol swab	1ea	Alcohol swab
2ea	Adhesive non-woven dressing with open cut in "U"shape	4ea	Gauze swab 7×7cm
4ea	Adhesive non-woven dressing for stabilizing the A-V lines	2ea	Surgical dressing in size 5×8cm
		2ea	Plastic band 40cm

#### Diabetes care



#### Safety insulin syringe







Safety insulin syringes with integrated needle. The innovative safety insulin syringe offers comfort and easy insulin therapy injection. And help eliminate accidental needle stick injury.

#### Features

- ▶ No gap design guaranteed the complete injection without dead space.
- ▶ Premium needle quality and silicon-coated to ease skin penetration and minimize pain.
- Convenient to use, safety mechanism will be activated when the intended dose has been
- ▶ Pull back the plunger after injection, the locking ring secures the needle in the barrel of the syringe preventing the reuse of syringes.
- ► Special design can prevent the needle back into the barrel occasionally.
- ► Visual confirmation of safety activation.
- ► Various range of syringes to meet the clinical requirements.



### Insulin syringe



#### Features

- ▶ Premium needle quality and silicon-coated to ease skin penetration and minimize pain.
- ► High transparency barrel with clearly readable graduated scale for better dose precision.
- ▶ No-dead-space.



Code	Unit	Size	Needle size
IS03	V U-100	0.3ml	25-33G×5/16"-1/2"
IS05	U-40/U-1000	0.5ml	25-33G×5/16"-1/2"
IS10	U-40/U-100	1ml	25-33G×5/16"-1/2"

Needle sizes according to customer's requests

#### Diabetes care



#### Safety insulin syringe







#### Features

- ▶ The innovative engineered safety mechanism to potentially reduce the risk of needle stick injury.
- ▶ By pull the safety cartridge to withdraw the needle into cartridge after injection.
- Easy to lock by rotating the safety cartridge.
- Low dead space minimizes medication waste.
- ► Triple sharpened needle guaranteed the sharpness.
- Clear graduation to guarantee the accuracy in small dose injection.



#### Using step

- Peel open the package and take out the syringe, ready for use.
  Extract the medication and pull up the safety cartridge to protect the needle, keep them at the surgical tray and ready for
- Pull back the safety cartridge and inject the drug medication by pushing the plunger.
  After the injection, pull up the safety cartridge again to cover the needle tip to avoid the injury risk.
- Rotating the safety cartridge until final locked with an audible "click" which activating the ultimate safety mechanism.
- ▶ Dispose the syringe in trash bin.

Code		Size	Needle size	
Type 1 /	Type 2			
ISS1-03	ISS2-03	0.3ml	29-33G×5/16"-1/2"	
ISS1-05	ISS2-05	0.5ml	29-33G×5/16"-1/2"	
ISS1-10	ISS2-10	1ml	28-33G×5/16"-1/2"	

#### Diabetes care



### Safety insulin pen needle









#### **Features**

- ► Automatically achieve the function of protecting needle tip after use.
- ▶ Obvious observation of use status.
- ► Easy self-disable set.
- ► Easy to see the droplets.

### Insulin pen needle



Thin wall technology ensures optimal insulin flow.

Premium needle quality and silicon-coated to ease skin penetration

Cartridge tip is anti-friction treated for easy, fragment-free cartridge puncture, preventing insulin contamination. Compatible with all pen injectors.

- ▶ NovoPen 4<sup>®</sup>
- ► NovoPen 3<sup>®</sup>
- ► NovoPen Junior®
- ► NovoRapidFlexpen®
- ▶ BD:Puregon Pen®
- ► Sanofi: Aventis Lantus®
- ► Sanofi: Optipen Pro 1<sup>®</sup>
- ► Sanofi: TactiPen®
- ► Owen Mumford: AutoPen®
  - ► Owen Mumford: AutoPen 2®
  - ▶ Lily: HumaPen®

Extra short needle length 5mm, 38% shorter than 8mm needle to guarantee the sufficient of insulin delivery and greater injection comfort.

Code	Gauge	Size and length (mm)	Code	Gauge	Size and length (mm)
Safety/Non-safety			Safety/Non-safety	MS	
IPNS33 IPN33	33G	0.2 x4-12mm	IPNŠ30 IPN30	30G	0.3 x4-12mm
IPNS32 IPN32	32G	0.23 x4-12mm	IPNS29 IPN29	29G	0.33 x4-12mm
IPNS31 IPN31	31G	0.25 x4-12mm			

### Diabetes care



### **Sharp disposal container for** insulin pen needle









#### Features

- ▶ The lid notches enable the pen needle to be removed safely, to avoid needle stick injury.

  Convenient for people handheld from all angles.
- No harmful medical grade material.
- ► Anti-leak flaps prevent leakage.



### **Alcohol swab**

Contain 70% isopropyl alcohol for optimum antibacterial action.



- Are thicker and softer than many other swabs, helping ensure reliable site preparation for safe, hygienic insulin injections. They are also appropriate for other household
- ▶ Individually foil wrapped in four layers.
- ▶ Provide an airtight seal that prevents leakage and drying.

